



# TECHNICAL SHEET

**Ikarios®**


Suture ligament anchor



## Suture ligament anchor

**Ikarios** is a non-absorbable suture intended for soft tissue, tendon and ligament reinsertion/repair.

IKARIOS is made from Ultra High Molecular Weight Polyethylene (UHMWPE) and is crimped with stainless steel needles, which are available in various sizes.



- + Implantable Medical Device in permanent contact with bone and soft tissues class IIb.
- + Ultra High Molecular Weight Polyethylene (UHMWPE) crimped with stainless steel needles.
- + Gamma ray sterilization.
- + Available in a variety of sizes with different needles.
- + Packaging by unit, recommended storage in a dry place protected from light. Period of validity of 60 months for sterile ones.
- + Informations in this brochure are intended to provide information on the product. Always refer to the packaging, labels and / or instructions for use before using this product. Document reserved for the use of health professionals.
- + This product is not available for sale in all countries.

CE 0459 - LNE GMED




### Product references

1266001	3 Metric (#2-0) - Length 45 cm 1 round needle
1266002	5 Metric (#2) - Length 93 cm 2 needles : 1 round et 1 triangular
1266003	7 Metric (#5) - Length 93 cm 2 triangular needles

Each reference is dedicated to a specific indication :

#### 1266001

For hand & foot

3/8  20mm  
Taper



#### 1266002

For knee, shoulder & foot

1/2  36.6mm  
Reverse cutting



1/2  25.9mm  
Taper



#### 1266003

For knee, shoulder, foot and hip

3/8  49 mm  
Sternotomy K



1/2  45 mm  
Sternotomy K

