



TECHNICAL SHEET


Foot implants

Compression and snap-off screws, varisation staple



Foot implants

Compression and snap-off screws, varisation staple



Compression screws for small bone osteosynthesis - mono and bi-cortical osteosynthesis, distal and proximal metatarsal osteotomies and hallux valgus osteotomies. Manufactured from titanium alloy (TA6V), available in different lengths with a 3mm diameter.

Snap off screws for small bone osteosynthesis – mono-cortical osteosynthesis and Weil osteotomies. Manufactured from titanium alloy (TA6V), available in different lengths with a 2mm diameter.

Varisation staples for realignment of small bone fragments – monocortical osteosynthesis of the phalange and phalangeal proximal or distal shortening osteotomies. Manufactured from stainless steel, available straight (90°) or oblique (26°).

- + Implantable Medical Device in permanent contact with bone and soft tissues > 30 days, class IIb.
- + Made of titanium TA6V for both screws, stainless steel for varisation staples.
- + Gamma ray sterilization - available non-sterile.
- + Available in a variety of sizes - see details in product references below.
- + Packaging by unit, recommended storage in a dry place protected from light. Period of validity of 60 months for sterile ones.
- + Informations in this brochure are intended to provide information on the product. Always refer to the packaging, labels and / or instructions for use before using this product. Document reserved for the use of health professionals.
- + This product is not available for sale in all countries.

CE 0459 - LNE GMED



Product references

Compression screw

110010	109010	Compression screw Ø 3 x 10 mm
110012	109012	Compression screw Ø 3 x 12 mm
110014	109014	Compression screw Ø 3 x 14 mm
110016	109016	Compression screw Ø 3 x 16 mm
110018	109018	Compression screw Ø 3 x 18 mm
110020	109020	Compression screw Ø 3 x 20 mm
110022	109022	Compression screw Ø 3 x 22 mm
110024	109024	Compression screw Ø 3 x 24 mm
110026	109026	Compression screw Ø 3 x 26 mm
110028	109028	Compression screw Ø 3 x 28 mm
110030	109030	Compression screw Ø 3 x 30 mm
110032	109032	Compression screw Ø 3 x 32 mm
110034	109034	Compression screw Ø 3 x 34 mm

Snap-off screw

110111	109111	Snap-off screw Ø 2 x 11 mm
110112	109112	Snap-off screw Ø 2 x 12 mm
110113	109113	Snap-off screw Ø 2 x 13 mm
110114	109114	Snap-off screw Ø 2 x 14 mm

Varisation staple

ST80	80	Straight varisation staple 8mm - 90°
ST81	81	Straight varisation staple 10mm - 90°
ST82	82	Oblique varisation staple 8mm - 26°
ST83	83	Oblique varisation staple 10mm - 26°

An instrument set is available to help the use of those implants :

