



TECHNICAL SHEET


Staples

Osteotomy



Staples

Osteotomy



The **osteotomy staples** are indicated for tibial valgus osteotomy (subtractive and additive).

Features an extra flat bridge (1.8 mm thickness) and rounded edges with no salient angle.

The staple driver securely holds the staple during insertion.

- + Implantable Medical Device in permanent contact with bone and soft tissues > 30 days, class IIb.
- + Made of chromium cobalt alloy
- + Gamma ray sterilization - available non-sterile.
- + Available in a variety of sizes, straight, with offset, semi-oblique or tripodial.
- + Packaging by unit, recommended storage in a dry place protected from light. Period of validity of 60 months for sterile ones.
- + Informations in this brochure are intended to provide information on the product. Always refer to the packaging, labels and / or instructions for use before using this product. Document reserved for the use of health professionals.
- + This product is not available for sale in all countries.

CE 0459 - LNE GMED

Product references

Straight osteotomy staples

(ST)20DP	H1 25 - L 20mm
(ST)26DP	H1 25 - L 26mm
(ST)36DP	H1 25 - L 36mm

Osteotomy staple with offset

(ST)73P	H1 36 - L 24 - D 5 - H2 25mm
(ST)74P	H1 40 - L 24 - D 10 - H2 25mm
(ST)75P	H1 44 - L 24 - D 14 - H2 25mm

Semi-oblique osteotomy staple

(ST)26SOP	H1 34 - L 26 - D 5 - H2 26mm
(ST)36SOP	H1 34 - L 34 - D 5 - H2 27mm

Tripodial osteotomy staple

26TSO	H1 34 - L 26 - D 5 - H2 27mm
-------	------------------------------

(ST) : available sterile or non sterile

Associated ancillary

0210-A	Staple driver
0300-A	Impactor
0400-A	Staple extractor with weight control
0410-A	Staple extractor

