



# TECHNICAL SHEET


## Interference screw

Metallic



# Interference screw

## Metallic



**METALLIC INTERFERENCE SCREW** is indicated for cruciate ligament surgery, tibial and femoral fixation in bone tunnels:

- of bone-tendon-bone graft for bone to bone interference screw
- of all types of grafts and of prosthetic Ligastic® ligaments for interference screws with atraumatic threads.

Their head has round edges and distal end is conical, ogive-shaped.

The screw is cannulated to allow the loading on a pin.

Their material is fully compatible with Magnetic Resonance Imaging (MRI).

An unique screwdriver fits all the sizes and make their implantation easier.

- + Implantable Medical Device in permanent contact with bone and soft tissues > 30 days, class IIb.
- + Made of titanium alloy (TA6V)
- + Gamma ray sterilization - available non-sterile.
- + Diameters 7, 8 or 9mm, length 18, 23 or 28mm
- + Packaging by unit, recommended storage in a dry place protected from light. Period of validity of 60 months for sterile ones.
- + Informations in this brochure are intended to provide information on the product. Always refer to the packaging, labels and / or instructions for use before using this product. Document reserved for the use of health professionals.
- + This product is not available for sale in all countries.

CE 0459 - LNE GMED

### Bone-to-bone

### Atraumatic



7 x 18 mm



7 x 18 mm



8 x 23 mm



8 x 23 mm



9 x 28 mm



9 x 28 mm



0010-C  
Metallic screws screwdriver

### Product references

Non sterile	Sterile	Description
<b>Bone-to-bone screws</b>		
7.18HEX	7.18HEX(S)	Interference screw Ø 7 x 18mm
8.18HEX	8.18HEX(S)	Interference screw Ø 8 x 18mm
9.18HEX	9.18HEX(S)	Interference screw Ø 9 x 18mm
7.23HEX	7.23HEX(S)	Interference screw Ø 7 x 23mm
8.23HEX	8.23HEX(S)	Interference screw Ø 8 x 23mm
9.23HEX	9.23HEX(S)	Interference screw Ø 9 x 23mm
7.28HEX	7.28HEX(S)	Interference screw Ø 7 x 28mm
8.28HEX	8.28HEX(S)	Interference screw Ø 8 x 28mm
9.28HEX	9.28HEX(S)	Interference screw Ø 9 x 28mm

<b>Atraumatic screws</b>		
0990718	1000718	Interference screw Ø 7 x 18mm
0990818	1000818	Interference screw Ø 8 x 18mm
0990918	1000918	Interference screw Ø 9 x 18mm
0990723	1000723	Interference screw Ø 7 x 23mm
0990823	1000823	Interference screw Ø 8 x 23mm
0990923	1000923	Interference screw Ø 9 x 23mm
0990728	1000728	Interference screw Ø 7 x 28mm
0990828	1000828	Interference screw Ø 8 x 28mm
0990928	1000928	Interference screw Ø 9 x 28mm

<b>Associated instrument</b>		
0010-C		Metallic screws screwdriver