



TECHNICAL SHEET


Staples

Ligament



Staples

Ligament



The **ligament staples** are intended for the fixation of autograft, allograft and prosthetic components.



Extra-flat bridge (1.8 mm thickness), rounded edges and without protruding angle.

They are designed with anti-back-out serrations with sharp bevelled points to prevent dislodging, depending on the reference.

The staple driver securely holds the staple during insertion.

- + Implantable Medical Device in permanent contact with bone and soft tissues > 30 days, class IIb.
- + Made of chromium cobalt
- + Gamma ray sterilization - available non-sterile.
- + Available in a variety of sizes, with or without spikes and with long or short legs.
- + Packaging by unit, recommended storage in a dry place protected from light. Period of validity of 60 months for sterile ones.
- + Informations in this brochure are intended to provide information on the product. Always refer to the packaging, labels and / or instructions for use before using this product. Document reserved for the use of health professionals.
- + This product is not available for sale in all countries.

CE 0459 - LNE GMED



Product references

Ligament staples with spikes

(ST)90P	Ligament staple 6 x 20mm
(ST)92P	Ligament staple 8 x 20mm
(ST)94P	Ligament staple 11 x 20mm
(ST)96P	Ligament staple 16 x 20mm

Ligament staples with spikes and short legs

(ST)90SP	Ligament staple 6 x 15mm
(ST)92SP	Ligament staple 8 x 15mm
(ST)94SP	Ligament staple 11 x 15mm
(ST)96SP	Ligament staple 16 x 15mm

Ligament staples without spike

(ST)91P	Ligament staple 6 x 20mm
(ST)93P	Ligament staple 8 x 20mm
(ST)95P	Ligament staple 11 x 20mm
(ST)97P	Ligament staple 16 x 20mm

Ligament staples without spike and short legs

(ST)91SP	Ligament staple 6 x 15mm
(ST)93SP	Ligament staple 8 x 15mm
(ST)95SP	Ligament staple 11 x 15mm
(ST)97SP	Ligament staple 16 x 15mm

Orengo ligament staple

(ST)ORENGO5	Oval staple 5 x 25mm
(ST)ORENGO10	Oval staple 10 x 25mm

(ST) : available sterile or non-sterile

Associated ancillaries

0210-A	Staple driver
0300-A	Impactor
0400-A	Universal staple extractor, control weight type
0410-A	Universal staple extractor
0251-A	Orengo staple driver